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Boston
Scientific
Advancing science for life™

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
VS

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Innova™ Vascular Self-Expanding Stent	(Boston Scientific)
Supera™ Peripheral Stent	(Abbott)
LifeStent™ Vascular Stent	(Bard)
EverFlex™ Self-Expanding Peripheral Stent	(Covidien/Medtronic)
S.M.A.R.T® Flex Vascular Stent and S.M.A.R.T. CONTROL® Vascular Stent	(Cordis/Cardinal)
Pulsar®-18	(Biotronik)
Compete® SE Vascular Stent	(Medtronic)

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
*Kaplan-Meier Estimate: Primary patency defined as core-lab assessed duplex ultrasound peak systolic velocity ratio (PSVR) ≤ 2.4 at 12 months in the absence of clinically-driven TLR or bypass of the target lesion.

**Log-rank p-value compares the entire K-M curves from time point zero to day 395 (full 1-year follow-up window)

1. EMINENT Trial: A global randomized controlled multi-center trial with 2:1 randomization of the Eluvia™ Drug-Eluting Stent against commercially-available Self-Expanding Bare Nitinol Stents, single-blind, superiority design; independent core lab adjudication. 12-Month Primary Patency rate of 83.2% in the Eluvia arm vs. 74.3% in the Bare-Metal Stenting arm (p -value = 0.0077).


EMINENT Clinical Trial 12-Month results presented by Professor Yann Goueffic, MD. VIVA 2021

Results from clinical studies are not predictive of results in other studies. Results in other studies may vary.

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1-YEAR PRIMARY SUSTAINED CLINICAL IMPROVEMENT ***

***In EMINENT, primary sustained clinical improvement was defined as an improvement (decrease) by at least 1 Rutherford category, without TLR.
EMINENT Clinical Trial 12-Month results presented by Professor Yann Goueffic, MD. VIVA 2021

Results from clinical studies are not predictive of results in other studies. Results in other studies may vary.

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As-treated. Major adverse events adjudicated by the Clinical Events Committee. P values from Chi-square test or two-sided Fisher's exact test.

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